Intramedullary Gliding Nails for children Intramed Gleitnagel für Kinder

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Description for use of Intramedullary Gliding Nails for Children:
For Shaft Fractures: Femoral, Tibial, Humeral and Ulna / Radius-Bone

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TREU INSTRUMENTE GmbH

- Traumatology – Orthopaedics

TREU manufactures since more than 4 decades trauma implants for bone healing, osteosyntheses products. Joint implants for hip surgery, Containers and Storage Cases and a Battery Driven Drill System are more major business fields.

To adapt to international standards and to achieve highest possible compatibility our standard osteosyntheses products are manufactured according to the guidelines of AO. Our Container and Storage Cases are oriented according to the world market leader.

Therefore, more than 8000 products which contribute worldwide, in more than 60 countries, to keep maintaining health and life style of millions of people stand in the center point of our entrepreneur ship.

In cooperation with leading surgeons we guarantee a continuous development and further development with the aim that with our products and services operation time can be reduced or healing process accelerated.

Since 1998 we are certified according to EU guidelines and manufacture our products in accordance with DIN - Specifications, ISO - Standards as well as US - and British Standards.

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Trocar with slotted sleeve for implantation of medullary gliding nails


After opening of the medullary canal with trocar, it happens quite often, that the point of entry cannot be found. By using the new trocar, the slotted sleeve remains in situ and avoids unnecessary searching or preparing of soft tissue. Gliding nails with bent tip can easily be introduced through the slot of the sleeve.

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

Opening of medullary canal with trocar and sleeve

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
Removal of trocar. Sleeve is held in position

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Introduction of nail through slot of sleeve

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The new Impactor for „Intramed Gliding Nails“ serves for quick introduction avoiding tiresome manipulating during fixation of implants. Owing to the lateral handle the danger of insurries to the surgeon in case of slipping-through of nails is practically impossible. A new clamping system serves for fixation of nails in the impactor without applying much power.

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9-1450 Pfriem
Perthes, 22,0 cm
9-1452 Pfriem
Perthes, 16,0 cm


9-1456 Pfriem
16,0 cm

1-5330 Extractions Forceps for Intramedullary
Gliding Nails
Extractionszange 22,0 cm für
Intramed-Schienen


7-388 Extraction Pliers for boring wire 130 mm
7-390 Extraction Pliers for boring wire 180 mm

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9-701 Slotted Hammer
Schlitzhammer

9-702 Bone Hammer / Weight approx: 350 g
Knochenhammer / Gewicht ca 350 gr.

7-434 Wire Cutter 220 mm TC
soft wire 3.0 mm
hard wire 2.5 mm

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
**Intramedullary Gliding Nail
PRODUCT NEWS**

**Intramed Gleitschiene
PRODUCT NEWS**

Intramedullary Gliding Nails

Diaphyseal fractures of long bones in children and adolescents in the growing age are mostly treated conservatively with plaster of Paris and external splintage. The borderline is there, where retention obviously can not be maintained. Quite often additional anaesthetics were necessary due to post-reduction. In spite of the often long therapeutical process the final result was not satisfactory neither for the child, nor for the parents and therapists. The treatment with plates and screws or with an external fixateur were reduced to exceptions.

Reasons for it:

- during growth 75% the periosteal healing is predominant
- with plates and screws a rigid anatomical reduction is given but followed by a long time of immobilization, extensive scar and reoperation due to metal removal.
- The treatment with external fixator for open fractures is advantageous and offers with good rigidity the possibility of early dynamization. However the axial positioning can not always be achieved. In addition regularly pin-care visits are necessary and the child always is confronted with the osteosynthesis material.

Demand

In spite of the „healing potential“ in childlike bone the treatment of bone in the growing skeleton the axial anatomical alignment must be looked at with priority.

Is there a reduction necessary then the fracture may not remain in an axial deviation or in a wrong rotational position. Independent from age the aim must be to achieve the optimum position.

Intramedullary osteosyntheses offers an alternative with

- little additional traumatization
- sufficient axial stability
- optimum healing due to stimulated micro movements
- low complication rate
- early mobilization
- excellent cosmetic result
- outpatient metal removal after 10 to 12 weeks and
- is indicated for patients, in which the conservative therapy has already started, but post reduction is predictable.


The use of about 10.000 intramedullary gliding nails in Germany has led to expand the indications also to adults.

This alternative to the treatment with intramedullary gliding nails results in:

- less invasive treatment
- shorter OR-time
- cost saving procedure
- less X-Ray exposure for patient and therapist
- outpatient metal removal after 6 to 9 months.


The rigidity of the elastic intramedullary gliding nail is adequate to the rigidity of the bone. Incoming forces (F_1), which lead to material breakage in rigid implants, are transformed and bend the implant. The elasticity of the implant and its ability of reversible deformation reset the intramedullary gliding nail into its origin position due to reset forces (R).

Pict. 1

If the intramedullary gliding nail (pict. 2) pre-bent (F_2), arises a new „Zero Point“, that means that the reset forces reset the gliding nail in its „new“ point of origin.


Pict. 2

Osteosyntheses with plates and screws (**pict. 3**) concentrate all biomechanical forces in a one point compression zone (A), so that bending forces on one side induces big tension forces on the opposite one (B).


With plates and screws:

Pict. 3

In case of a rigid intramedullary nail (**pict. 4**), which clamps inside the intramedullary canal, changing induced forces may lead to varus deviation of the bone due to implant deforming or may result in nail breakage (D).

Rigid Nail


Pict. 4

Pict. 5


The basic principle are two intramedullary gliding nails in opposite position which form a Y-shaped buttress (A). The induced forces on the implant are reduced at about 30%.

Pict. 5

Pict. 6

The intramedullary gliding nail may absorb incoming forces, which lead to implant breackage in rigid implants, due to the ability of the reversable deforming.


Pict. 6

Pict. 7


Axial bending forces (F) distract a gliding nail. Its reset-forces (R) correct the axial failure.

Pict. 7

Pict. 8


Frontal pushing forces distract both gliding nails. The appropriate reset-forces push them back into its origin position.

Pict. 8

Pict. 9


In oblique fractures the axial induced force (F) is split into a pushing- (S) and a compression-component (K). The intramedullary guiding avoids (S) and strengthens (K).

Pict. 9

Pict. 10

Rotational winding of the gliding nails one around each other built up tension to correct malpositioning.


Pict. 10

Pict. 11


Sagittal pushing forces reduce the contact-area within the intramedullary canal and provoke, due to induced tension, the reset of the fragments. Are there excentric sagittal forces it leads to ante- or recurvation which will be corrected by the elasticity of the implants.

Pict. 11

Pict. 12

Due to the fact that the oval between the two gliding nails cannot be enlarged, axial compression leads to re-inforcement pressure on the gliding nails against the endost.

Pict. 12

Technique:

Determination of the thickness size of the gliding nail

$$\text{Gliding nail size} = \frac{\text{Medullary canal diameter in mid-shaft}}{3}$$

$$\text{Forearm} = \frac{2 \times \text{Medullary canal diameter in mid-shaft}}{3}$$

Entry point:

The entry point for gliding nails must always be outside the joint capsule and care to be taken about the epiphyseal area. A region with little soft tissue coverage and appropriate intramedullary canal is advantageous.

Forearm

Always one gliding (**pict. 13**) nail is used in a conventional way starting from distal radius or from proximal ulna.

Pict. 13

Humerus

The ascending humeral gliding nails (**pict. 14**) can be performed either from both sides or both nails are inserted from radial (**pict. 15**) approach. In this case There is one incision but two drill holes.

Pict. 14

Pict. 15

Pict. 16
Femur fractures routinely are treated with the both side distal approach.


Pict. 16

Pict. 17
Distal humeral- and femoral fractures in some cases need an ascending treatment. There is one incision but two drill holes.

Pict. 17

Pict. 18
Tibia fractures are treated in a descending manner with two incisions at both sides of tuberositas tibiae.

Pict. 18


Pict. 19

OP- Technique

Skin incision is made from the planned entrance point in epiphyseal direction (**pict. 19**). A 2 to 3 cm long incision facilitates the following steps.


Pict. 20


Pict. 21


The entrance point is marked with the trokar which has been put into the sleeve to protect soft tissue. The marking is done in a 90° angulation (**pict. 20**) to the bone surface. The bone opening can alternatively done either with the trokar or with a drill (**pict. 21**).

Pict. 22

Pict. 22a

Pict. 23

The size of the trokar or drill must be at least 0,5 mm bigger in diameter than the indicated intramedullary nail. The bone perforation is performed in angulation less than 45° (**pict. 22**). After removal of the trokar or drill the tissue protection sleeve is held in position (**pict. 22a**). The gliding nail is then introduced through it. The slot of the sleeve (**pict. 23**) facilitates the introduction of nails with bent tips.


Implantation of intramedullary gliding nail

Is pre-bending of the gliding nail necessary? It is then indicated when the gliding nail in the bone entrance fragment shall early reach the opposite cortical bone, e.g. in those cases when the fragment is relatively short or if the implant sticks on the opposite cortical bone and cannot be moved further on into the medullary canal. Primarily that implant is introduced which leads to the highest reduction effect of the fracture. The gliding nail is introduced through the tissue sleeve into the intramedullary canal. The gliding nail then must be turned in the way that the tip shows in direction to the medullary canal and not to the cortex.

Pict. 24


With the T-handle with Jacobs Chuck the gliding nail is shortly fixed and introduced step by step. This facilitates the introduction and avoids unintended bending of the implant.

Pict. 24

Pict. 25

Is there little fragment contact at the fracture side, the tip of the gliding nail can be aligned while twisting the gliding nail.

Pict. 25


Tip:

Often it is useful, that prior to the transfer of the first gliding nail from one side of the fracture to the other one that the second gliding nail has been inserted and being pushed onto the borderline of the fracture. There are now two starting points given which can be used for reduction in detail.

Pict. 26


To introduce the second gliding nail in the opposite fragment easily, the first gliding nail can be twisted to correct fragment position.

Pict. 26

Pict. 27


The implant is fixed in the metaphyseal spongiosa of the counter fragment with some soft impacts with the hammer.

Pict. 27

Pict. 28

Slight misalignments in axial position can be corrected by restricted twisting manoeuvres of the gliding nail.

Pict. 28


Tip:

A possible distraction at the fracture line which may occur due to the implantation can be corrected through axial compression. Afterwards the implants are impacted. Finally the gliding nail will be cut subcutaneously with the wire cutter. A protective cap is put onto the sharp end of the gliding nail to prevent soft tissue irritation and nail perforation.

Pict. 29

In oblique fractures that gliding nail is easier to implant which nail tip is 90° to the fracture area.

Pict. 29

Pict. 30

In spiral fractures primarily the implantation at the side with the long cortical bone is recommended.

Pict. 30

Pict. 31

Abb. 32

Implant removal

Implant removal can be done in outpatient treatment. The end of the gliding nail is exposed in the conventional way. The protection cap is removed with a forceps (**pict. 31**). The implant is grasped with an adjustable extraction forceps (**pict. 32**), with a longitudinal groove to cover the implant.

Implantation Problems

Pict. 33

The assymetrical implantation leads to two tension bows and therefore can induce malposition.

Pict. 33

Pict. 34

Too short or too thin used implants are not able to fulfil the demands of the elastic-stable treatment with intramedullary gliding nails.

Pict. 34

Pict.35

Twisted implants one around eachother avoid the adequate tension fixation and the necessary elasticity.

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